PJX Series Balances
Instruction Manual

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Table of Contents

1. INTRODUCTION	EN-3
1.1 Description	EN-3
1.2 Features	EN-3
1.3 Definition of Signal Warnings and Symbols	EN-3
1.4 Safety Precautions	EN-3
2. INSTALLATION	EN-4
2.1 Unpacking	EN-4
2.2 Select the Location	EN-4
2.3 Leveling	EN-4
2.4 Connecting Power and Acclimating the Balance	EN-4
2.5 Connecting the Interface	EN-5
2.6 Initial Calibration	EN-5
3. OPERATION	EN-6
3.1 Overview of Display, Home Screen	EN-6
3.2 Principal Functions and Main Menu	EN-7
3.3 Overview of Parts and Features – Draft Shield Models	EN-7
3.4 Overview of Parts and Features – Non-Draft Shield Models	EN-7
4. APPLICATIONS	EN-8
4.1 Weighing	EN-8
4.2 Parts Counting	EN-8
4.3 Percent Weighing	EN-10
4.4 Density Determination	EN-12
4.4.1 Measuring the Density of a Sinking Solid Using Water	EN-16
4.4.2 Measuring the Density of a Floating Solid Using Water	EN-17
4.5 Additional Features	EN-18
5. MENU SETTINGS	EN-19
5.1 Menu Navigation	EN-19
5.1.1 Changing Settings	EN-19
5.2 Calibration	EN-19
5.2.1 Calibration Sub-menu (InCal models)	EN-19
5.2.2 Internal Calibration (not applicable to ExCal models)	EN-19
5.2.3 Cal Adjust (not applicable to ExCal models)	EN-20
5.2.4 Span Calibration	EN-20
5.2.5 Linearity Calibration	EN-21
5.3 Balance Setup	EN-23
5.3.1 Filter Level	EN-23
5.3.2 AZT (Auto Zero Tracking)	EN-23
5.3.3 Auto Tare	EN-23
5.3.4 Graduations	EN-24
5.3.5 Date Format	EN-24
5.3.6 Date Setup	EN-24
5.3.7 Time Format	EN-24
5.3.8 Time Setup	EN-25
5.3.9 Brightness	EN-25
5.3.10 Auto Dim	EN-25
5.3.11 Approved Mode	EN-25
5.4 Weighing Units	EN-26
5.5 RS232 Interface Setup	EN-27
5.5.1 Baud Rate	EN-27
5.5.2 Transmission	EN-28
5.5.3 Handshake	EN-28
5.6 Print Settings	EN-28
5.6.1 Stable Only	EN-28
5.6.2 Numeric Only	EN-29
5.6.3 Single Header	EN-29
5.6.4 Auto Print	EN-29
5.6.5 Header	EN-30

5.6.6	Date and Time.....	EN-30
5.6.7	Balance ID.....	EN-30
5.6.8	Balance Name.....	EN-30
5.6.9	User Name.....	EN-30
5.6.10	Project Name.....	EN-30
5.6.11	Application Name.....	EN-30
5.6.12	Result.....	EN-30
5.6.13	Gross.....	EN-30
5.6.14	Net.....	EN-31
5.6.15	Tare.....	EN-31
5.6.16	Line Feed.....	EN-31
5.7	GLP.....	EN-31
5.7.1	Header.....	EN-31
5.7.2	Balance Name.....	EN-32
5.7.3	User Name.....	EN-32
5.7.4	Project Name.....	EN-32
5.8	Factory Reset.....	EN-32
5.9	Lockout.....	EN-32
6.	LEGAL FOR TRADE (LFT)	EN-33
6.1	Settings.....	EN-33
6.2	Verification.....	EN-33
6.3	Sealing.....	EN-33
7.	PRINTING	EN-34
7.1	Connecting, configuring and Testing the Printer / Computer Interface.....	EN-34
7.2	Output format.....	EN-35
7.3	Printout Examples.....	EN-35
8.	MAINTENANCE	EN-37
8.1	Calibration.....	EN-37
8.2	Cleaning.....	EN-37
8.3	Troubleshooting.....	EN-37
8.4	Service Information.....	EN-37
9.	TECHNICAL DATA.....	EN-38
9.1	Specifications.....	EN-38
9.2	Drawings and Dimensions.....	EN-44
9.3	Accessories.....	EN-44
9.4	Communication.....	EN-45
9.4.1	Interface Commands.....	EN-45
9.4.2	RS232 (DB9) Pin Connections.....	EN-46
10.	SOFTWARE UPDATES.....	EN-46
11.	COMPLIANCE	EN-47

1. INTRODUCTION

1.1 Description

The PJX balance is a precision weighing instrument that will provide you with years of service if properly cared for. PJX balances are available in capacities from 40 grams to 6200 grams.

1.2 Features

Operation Controls: backlit LCD display, with 4 weighing applications and many features.

1.3 Definition of Signal Warnings and Symbols

Safety notes are marked with signal words and warning symbols. These show safety issues and warnings. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results.

- WARNING** For a hazardous situation with medium risk, possibly resulting in injuries or death if not avoided.
- CAUTION** For a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or injuries if not avoided.
- Attention Note** For important information about the product
For useful information about the product

Warning Symbols

General Hazard

Electrical Shock Hazard

Alternating current

Direct current

1.4 Safety Precautions

CAUTION: Read all safety warnings before installing, making connections, or servicing this equipment. Failure to comply with these warnings could result in personal injury and/or property damage. Retain all instructions for future reference.

- Verify that the AC adapter's input voltage range and plug type are compatible with the local AC mains power supply.
- Make sure that the power cord does not pose a potential obstacle or tripping hazard.
- Do not position the balance such that it is difficult to reach the power connection.
- The balance is for indoor use only. Do not operate the equipment in hazardous or unstable environments.
- Operate the equipment only under ambient conditions specified in these instructions.
- Do not drop loads on the pan.
- Use the balance only in dry locations.
- Disconnect the equipment from the power supply when cleaning.
- Use only approved accessories and peripherals.
- Service should only be performed by authorized personnel.

2. INSTALLATION

2.1 Unpacking




Carefully remove your PJX balance and each of its components from the package. The included components vary depending on the balance model. Save the packaging to ensure safe storage and transport. Please read the manual completely before installing and using the PJX balance to avoid incorrect operation.

Components included:

- Balance
- Power adapter + Attaching plug
- Stainless steel pan
- Pan support (for 0.1 g / 0.01 g model only)
- Warranty card

2.2 Select the Location

Avoid heat sources, rapid temperature changes, air current or excessive vibrations. Allow sufficient space.


2.3 Leveling

Be sure the balance is level before it is used or after its location is changed.

The PJX balance has a level bubble in a small round window beside the display.

To level the balance, adjust the 4 Leveling Feet until the bubble is centered in the circle.

Please refer to the right figure for leveling.

2.4 Connecting Power and Acclimatising the Balance

Connect the DC output connector to the power receptacle on the rear of the balance. Then connect the AC adapter plug to a suitable electrical outlet.

Acclimatising

It is suggested that the balance should not be used until it has been connected to power and acclimatised to the environment for a certain period of time. In the case of a balance with the precision above 0.1 mg, the acclimatisation time should be 1.5 hours; in the case of balance with the precision of 0.01 mg, the acclimatisation time should be more than 4 hours.

2.5 Connecting the Interface

The PJX balance has two connectivities, RS232 port and Mini DIN.

Use the RS-232 port to connect either to a computer or a printer with a standard (straight-through) serial cable.

Mini DIN is used to connect with another accessory, as an alternative option for RS232.

Interface connections on the rear of the balance.

Mini DIN

RS232

RS232: Used to connect to PC or Printer

Note: See the Printing section for Connecting, Configuring and Testing the Printer / Computer Interface.

2.6 Initial Calibration

When the PJX balance is first installed, or when it is moved to another location, it must be calibrated to ensure accurate weighing results. PJX balances are classified into two categories, InCal models and ExCal models. InCal models have a built-in calibration mechanism which can calibrate the balance automatically and does not require the use of external calibration masses. If preferred, InCal models can also be manually calibrated with external masses. ExCal models are calibrated with external masses. Make sure to have the appropriate calibration masses available before beginning calibration.

3. OPERATION

3.1 Overview of Controls and Display

CONTROLS

CONTROL FUNCTIONS

Button	 Zero Yes	 Print Unit No	 Function Mode Back	 Tare Menu-Cal Exit
Primary Function (Short Press) 	On / Zero <ul style="list-style-type: none"> If the balance is Off, turns on the balance. If balance is On, sets zero. 	Print <ul style="list-style-type: none"> Sends the current displayed value to the serial interface. 	Function <ul style="list-style-type: none"> Operation is dependent on the application mode. 	Tare <ul style="list-style-type: none"> Performs tare operation.
Secondary Function (Long Press) 	Off <ul style="list-style-type: none"> Zeroing current value. 	Unit <ul style="list-style-type: none"> Changes weighing units. 	Mode <ul style="list-style-type: none"> Changes application mode. 	Menu-Cal <ul style="list-style-type: none"> Enters the main menu. Calibration is the first sub-menu. Views the preset Tare value.
Menu Function (Short Press) 	Yes <ul style="list-style-type: none"> Accepts the current (blinking) setting on the display. 	No <ul style="list-style-type: none"> Rejects the current (blinking) setting on the display. Increments a value being entered. 	Back <ul style="list-style-type: none"> Reverts back to previous menu item. Decrements a value being entered. 	Exit <ul style="list-style-type: none"> Immediately exits the sub-menu. Aborts a calibration in progress.


MAIN APPLICATION SCREEN


3.2 Principal Functions and Main Menu

- Weighing:** First press **Zero** to set the display to zero. Place an object on the pan. The display indicates the gross weight.
- Taring:** With no load on the pan, press **Zero** to set the display to zero. Place an empty container on the pan and press **Tare**. Add material to the container and its net weight is displayed. After the container and the objects are removed, the load will be displayed as a negative number. Press **Tare** to clear.
- Zero:** Press **Zero** to zero the balance.

3.3 Overview of Parts and Features – Draft Shield Models

3.4 Overview of Parts and Features – Non-Draft Shield Models

4. APPLICATIONS

The PJX balance can be operated in 4 application modes by long pressing the **Function / Mode** button.

4.1 Weighing

Note: Before using any application, be sure the balance has been leveled and calibrated.

Use this application to determine the weight of items in the selected unit of measure.

Weighing

<ol style="list-style-type: none"> 1. Press Tare or Zero if necessary to begin. 2. Press and hold the Function / Mode button to select WEIGH (this application is the default). 	
<ol style="list-style-type: none"> 3. Place objects on the pan to display the weight. Once the reading is stable, the * will appear. 4. The resulting value is displayed in the active unit of measure. 	

Item Settings

To view or adjust the current settings.


- **Weighing Units:** Change the displayed unit. See Section 5.4 for the detailed processes.
- **Filter Level:** Change Filtering level. See Section 5.3.1 for more information.
- **GLP Data:** See Section 5.7 for more information.
- **Print Settings:** Change printing settings. See Section 7 for more information.



4.2 Parts Counting

Note: Before using any application, be sure the balance has been leveled and calibrated. The minimum piece weight should be no less than 0.1d.

Use this application to count samples of uniform weight.

Parts Counting

<ol style="list-style-type: none"> 1. Press Tare or Zero if necessary to begin. 2. Press and hold the Function / Mode button until Count appears. 	
---	--

<p>3. After confirmation by pressing Yes, the message CLr.APW will appear on the screen.</p>	
<p>4. Press Yes, and the message Pwt will display with the numeral 15 (default) flashing. The user can press No or Back to increase or decrease the value. For instance, to increase the value to 15, please press Yes. Then, Pwt and 15 will flash simultaneously.</p>	
<p>5. Place 15 samples on the pan. Press the Function / Mode button so that the weight of the 15 samples is used to establish the average piece weight (APW). The display will show 15 pieces.</p>	
<p>6. Remove the 15 samples from the pan and then place additional samples on the pan. The corresponding number of pieces will display on the screen.</p>	
<p>7. To view the total weight or the number of pieces of the objects, press the Function / Mode button.</p>	

Item Settings

To view or adjust the current settings.

<p>Sample size: The sample size ranges from 1 to 100. The default value is 10.</p> <p>Note: If the APW of the last parts counting operation needs to be kept, the user can press No when the display shows the message <code>CLr.APW</code> (clear the average piece weight. Place additional objects on the pan, and the corresponding number of pieces will display.</p>	
<p>APW Optimization: Improving counting accuracy by re-calculating the piece weight automatically as parts are added.</p> <p>APW Optimization occurs only when the number of pieces added to the pan is between one and three times the number already on the pan.</p>	
<p>Print Settings: Changing printing setup. See Section 7 for more information.</p>	


4.3 Percent Weighing




Note: Before using any application, be sure the balance has been leveled and calibrated.

Use Percent Weighing to display the weight of a test object as a percentage of a pre-established reference sample.

Note: The minimum reference weight should be no less than 0.1d.


Percent Weighing

<p>1. Press and hold the Function / Mode button until <code>PERCENT</code> appears.</p>	
<p>2. After confirmation by pressing Yes, the message <code>CLr.rEF</code> (clear the reference) will appear on the screen.</p>	

<p>3. Press Yes, and <i>PuL.rEF</i> (put the reference weight) will display.</p>	
<p>4. Place the reference sample on the pan to display the weight. When the reading is stable, the * appears. 5. Press the Function / Mode button so that the weight of the reference sample is stored in memory. The display will show <i>100%</i>.</p>	
<p>6. Remove the reference sample, and place the test object on the pan. The ratio of the test object to the reference sample weight is displayed as a percentage.</p>	
<p>7. To view the reference sample weight or the test object weight, press the Function / Mode button.</p>	
<p>8. To establish a new reference sample weight, press and hold the Function / Mode button and repeat the steps described above.</p>	

Item Settings

Note: If the reference weight of last Percent Weighing operation needs to be kept, press **No** when the message **CLr.rEF** (Clear reference) displays.

Printing Setup:

Changing printing setup. See Section 7 for more information.


4.4 Density Determination

Note: Before using any application, be sure the balance has been leveled and calibrated.

Use this application to determine an object's density.

A **Density Determination Kit, Part Number 80253384**, is designed to be used with PJX series balances. Illustrations in this procedure refer to the density kit, however, you may use whatever lab apparatus that will suit the requirements for density measurements. A built in reference density table for water at the temperature of 20°C is included in the balance software. Review this entire section before attempting density measurements.

When making density measurements, the material should weigh at least 10.0 mg on an analytical balance and 100 mg on a precision balance.


Balance Preparation with Ohaus Density Kit (Optional)

Allow the balance to warm up sufficiently before making measurements.



Open either the left or right side door of the balance and remove the Pan as shown. Insert the Bracket into the balance where the Pan was removed.

The Equalizing Washer is not used.

- 1. Glass beaker
- 2. Bracket
- 3. Forceps
- 4. Platforms
- 5. Holder for floating solids
- 6. Holder for non floating solids
- 7. Precision thermometer with holder
- 8. Sinkers 10ml (optional equipment)

Place the Support into position over the bracket making sure the Support does not make contact with the Bracket as shown in illustration.

Install beaker on support as shown.
NOTE: Beaker and thermometer are not supplied as part of the density kit.

- The density Q is the quotient of the mass m and the volume V.

$$Q = \frac{m}{V}$$

Density determinations are performed by using Archimedes' principle. This principle states that every solid body immersed in a fluid loses weight by an amount equal to that of the fluid it displaces. The density table for water is included in the Discovery balance software.


The density of a solid is determined with the aid of a liquid whose density, Q₀, is known (water is used as an auxiliary liquid). The solid is weighed in air (A) and then in the auxiliary liquid (B). The density Q can be calculated from the two weighings as follows:

$$Q = \frac{A}{A - B} \cdot Q_0$$


The balance allows direct determination of the buoyancy P (P = A - B) and consequently the above formula can be simplified:

$$Q = \frac{A}{P} \cdot Q_0$$


- Q = Density of the solid
- A = Weight of the solid in air
- B = Weight of the solid in the auxiliary liquid
- Q₀ = Density of the auxiliary liquid at a given temperature (this value depends on the temperature). The density table for water is included in Discovery balances.
- P = Buoyancy of the solid in the auxiliary liquid (corresponds to A-B).

Sample Weighing in Air

Sample Weighing in Liquid

Buoyancy Sample Weighing

Place the solid in the Weighing Pan on the Weigh Below Hook in the liquid as shown. Ensure that there are no air bubbles on the solid to be weighed. Close the draft shield doors and weigh the solid (buoyancy P). The display indicates the density in grams/cc.

Solid Density Determinations for items Less Dense Than Water

For density determination of solids with a density less than 1 g/CM³, the bottom of the Weigh Below Hook for solids must be used as it holds the solid body below the surface of the auxiliary liquid. If the buoyancy of the solid is greater than the weight of the Weigh Below Hook, the Weigh Below Hook must be weighted by placing an additional mass on the submerged part of the Weigh Below Hook as shown. Weigh the sample in air first as explained in the previous procedure. After loading the additional mass, tare the balance and start the weighing again. Wait until the balance has reached stability and note the displayed weight P (buoyancy of the solid).

Improving the Accuracy of the Result of Solid Density

The following tips should help you improve the accuracy of the results in the density determination of solids.

Temperature

Solids are generally so insensitive to temperature fluctuations that the corresponding density changes are of no consequence. However, as work is performed with an auxiliary liquid in the density determination of solids, their temperature must be taken into account as the temperature has a greater effect with liquids and causes density changes in the order of magnitude 0.1 to 1% per °C. This effect is already apparent in the third decimal place of the result.

To obtain accurate results, we recommend that you always take the temperature of the auxiliary liquid into account on all density determinations.

Air Buoyancy

1 CM³ of air weighs approximately 1.2 mg (depending on the physical condition). As a consequence, in the weighing in air, each solid experiences buoyancy of this magnitude (the so-called "air buoyancy") per cm³ of its volume.

However, the air buoyancy must be taken into account only when a result is required with an accuracy of 3 to 4 decimal places. To correct for this, the air buoyancy (0.0012 g per cm³ volume of the body) is added to the calculated result:

$\text{Calculated density} + 0.0012 \text{ g/cm}^3 \text{ air buoyancy} = \text{effective density}$

Surface tension of the auxiliary liquid

Adhesion of the liquid to the Weigh Below Hook causes an apparent weight increase of up 3 mg.

As the Weigh Below Hook is immersed in the auxiliary liquid in both weighings of the solid (in air and in the auxiliary liquid), the influence of the apparent weight increase can be neglected because the balance is tared before every measurement.

To reduce the effect of air bubbles and to ensure the greatest possible accuracy, use a few drops of a wetting agent (not supplied) and add them to the auxiliary liquid.

Liquid Density Determinations



The density of a liquid can be made using a sinker of known volume. The sinker (P/N: 83034024) is weighed in air and then in the liquid whose density is to be determined, The density, Q, can be determined from the two weighings as follows:

$Q = \frac{A - B}{V}$

- Q = Density of the liquid
- A = Weight of the sinker in air
- B = Weight of the sinker in liquid
- V = Volume of the sinker
- P = Buoyancy of the sinker in the liquid (P = A-B)

In DENSITY SETUP, set Mode to Liquid Density and enter sinker volume in cc's.

After weighing the sinker in air and then weighing the sinker immersed in liquid, the balance calculates the density of the liquid and is displayed in grams/cc. See illustrations below for placement of the sinker. When the sinker is immersed in the liquid, it must not come into contact with the bottom of the beaker.




Operation Method


Two types of density determination can be made:

1. Solid more dense than the water
2. Solid less dense than the water

Note: It is required that the auxiliary liquid should be water and the water temperature should be kept at 20°C.

4.4.1 Measuring the Density of a Sinking Solid Using Water

<p>1. Press and hold the Function / Mode button until the dENs itY appears on the screen.</p>	
<p>2. After confirmation by pressing Yes, the message Air Luft will display on the screen.</p>	
<p>3. Weigh the sample in air using the balance and the density determination kit. When the * (symbol of stability) appears, press the Function / Mode button to confirm the weight of the sample in the air.</p>	
<p>4. Weigh the sample submerged in the liquid using the balance and the density determination kit. Note: Push the sample down into the liquid until it is fully submerged.</p>	
<p>5. Press the Function / Mode button to get the density of the sample.</p>	

<p>6. After completion of the test, press the Function / Mode button to test a new sample.</p>	 <p>The image shows the Pioneer balance display with the number 1.9938. Below the display are four buttons: 'Zero' (Yes), 'Print Unit' (No), 'Function Mode' (Back), and 'Tare Menu-Cal' (Exit). A hand icon is shown pressing the 'Function Mode' button.</p>
---	--

4.4.2 Measuring the Density of a Floating Solid Using Water

<p>Long press the Function / Mode button until the $dEN5.1tY$ appears on the screen.</p>	 <p>The image shows the Pioneer balance display with the text dEN5.1tY. Below the display are four buttons: 'Zero' (Yes), 'Print Unit' (No), 'Function Mode' (Back), and 'Tare Menu-Cal' (Exit). A hand icon is shown long-pressing the 'Function Mode' button.</p>
--	---


<p>In determining density with the balance, the balance setup and density determination procedures are basically the same for a floating solid and a non-floating solid except for the necessary holder (as shown in the figure) to be used for density determination.</p> <p>After completion of the test, press the Function / Mode button to test a new sample.</p> <p>Note: Push the sample down into the liquid until it is fully submerged.</p>	 <p>A diagram of a cylindrical holder with a flat bottom and a narrow neck, used for non-floating solids.</p>	 <p>A diagram of a cylindrical holder with a flat bottom and a narrow neck, containing a liquid and a solid sample floating on the surface.</p>
	<p>Holder for non-floating solids</p>	<p>Holder for floating solids</p>

4.5 Additional Features

Weigh Below

Note: Ensure the balance has been leveled and calibrated.

The PJX balance is equipped with a weigh below hook for weighing below the balance (as shown in the figure below).

Before turning the balance over, remove the pan and draft shield elements (if present) to prevent damage. Do not place the balance on the pan support cone or load cell pins.

To use this feature, remove power from the balance, then remove the protective cover for the weigh below opening.

Power on the balance, and then use a string or wire to attach items to be weighed.

Weigh below protective cover

Weigh below hook

5. MENU SETTINGS

5.1 Menu Navigation

Calibration	Setup	Unit	RS232	Print	GLP	Reset	Lock
InCal	Filter Level	Gram	Baud Rate	Stable Only	Header 1	Reset All	Calibration
Cal Adjust	AZT	Kilogram	Parity	Numeric Only	Header 2		Setup
Span Cal	Auto Tare	Milligram	Handshake	Single Header	Header 3		RS232
Linearity Cal	Graduations	Carat		Auto Print	Balance Name		Print
	Date Format	Pound		Header	User Name		GLP
	Date Setting	Ounce		Date and Time	Project Name		Reset
	Time Format	Ounce Troy		Balance ID			
	Time Setting	Pennyweight		Balance Name			
	Brightness	Grain		User Name			
	Auto Dim			Project Name			
	LFT			Application Name			
				Result			
			Gross Weight				
			Net Weight				
			Tare Weight				
			Line Feed				

5.1.1 Changing Settings

To change a menu setting, navigate to that setting using the following steps:

Enter the Menu

Long press the Menu button to enter the **Menu**.

Select the Sub-Menu

Press **No** to step between the sub-menus, and press **Yes** to enter the sub-menu.

Select the Menu Item

Press **No** to step through the Menu Items, and press **Yes** to enter the displayed Menu Item.

5.2 Calibration

PJX balances offer a choice of three calibration methods: Internal Calibration (for InCal models only), Span calibration and Linearity Calibration.

Attention: Do not disturb the balance during any calibration.

5.2.1 Calibration Sub-menu (InCal models)

Note: ExCal models only have Span Calibration and Linearity Calibration.

5.2.2 Internal Calibration (not applicable to ExCal models)

Calibration is accomplished with the internal calibration mass. Internal Calibration can be performed at any time, provided the balance has warmed up to operating temperature and is level.

With the Balance turned ON and no load on the pan, press the **Tare / Menu-Cal** button and select **InCal** to initiate the internal calibration.

The screen shows the status, and then press any button to return to the current application after calibration.

5.2.3 Cal Adjust (not applicable to ExCal models)

Use this calibration method to fine tune the effect of the Internal Calibration.

Calibration Adjust may be used to adjust the result of the Internal Calibration by ± 100 divisions.

Note: Before making a calibration adjustment, perform an Internal Calibration. To verify whether an adjustment is needed, place a test mass equal to the **span calibration value** on the pan and note the difference (in divisions) between the nominal mass value and the actual balance reading. If the difference is within +/- division, calibration adjustment is not required. If the difference exceeds +/-1 division, calibration adjustment is recommended.

Example:

Expected weight reading:	200.000g (Test mass value)
Actual weight reading:	200.014g
Difference in grams:	- 0.014g
Difference in divisions:	- 14 (InCal Adjust value)

To perform a Calibration Adjustment, select InCal Adjustment from the list of Calibration Menu; enter the value (positive or negative divisions) to match the difference noted earlier in the procedure.



Recalibrate using Internal Calibration. After calibration, place the test mass on the pan and verify that the mass value now matches the displayed value. If not, repeat the procedure until Internal Calibration reading agrees with the test mass.

5.2.4 Span Calibration

Span calibration uses two calibration points, one at **zero load** and the other at **specified full load** (span). For detailed calibration mass information please refer to the specification tables in the "Span Calibration Points", SPECIFICATIONS, Section 9.

With the balance turned On and no load on the pan, Span Calibration can be performed. The best accuracy is achieved using the mass closest to the full span value.

Steps for span calibration


<p>1. Press and hold the Tare / Menu-Cal button, and the <i>Calibration Menu</i> will display.</p>	
<p>2. Press Yes to enter the Calibration Menu.</p>	



<p>3. To change the calibration mode, press No until SPAN (span calibration) is displayed.</p>	
<p>4. The Calibration mass value will be shown in the screen. After the display shows 200.0000 g, place weight(s) of 200 g on the pan for calibration. To change to the calibration point of half full capacity (e.g. 100 g), press the Function / Mode button. After the screen shows "Place weights" and "100.000 g", place weight(s) of 100 g on the pan for calibration.</p>	
<p>5. Once the span calibration is completed successfully, CALdone will display. Press any button to return to the previous screen.</p>	
<p>6. Remove the weight(s), and the reading will be set to zero.</p>	

5.2.5 Linearity Calibration

Linearity calibration uses three calibration points, one at zero load and the others at specified loads. With no load on the balance, press Linearity Calibration to begin the process. The balance captures the zero point, and then prompts for the next weight. Continue to follow the instructions on the display until the calibration is completed.

Steps for linearity calibration

<p>1. Long press the Tare / Menu-Cal button, and the Calibration Menu will display.</p>	
--	--

<p>2. Press Yes to enter the Calibration Menu.</p>	 <p>PIONEER™</p> <p>CAL</p> <p>Zero  Yes </p> <p>Print Unit No</p> <p>Function Mode Back</p> <p>Tare Menu-Cal Exit</p>
<p>3. To change the calibration mode, press No until LINEAR (linearity calibration) is displayed.</p>	 <p>PIONEER™</p> <p>LINEAR</p> <p>Zero  Yes</p> <p>Print Unit No </p> <p>Function Mode Back</p> <p>Tare Menu-Cal Exit</p>
<p>4. The calibration mass value will be shown in the display. After the display shows 100.0000 g, place weight(s) of 100 g on the pan for calibration.</p>	 <p>PIONEER™</p> <p>100.0000 g</p> <p>Zero  Yes</p> <p>Print Unit No</p> <p>Function Mode Back</p> <p>Tare Menu-Cal Exit</p>
<p>5. Remove the weight(s) of 100 g from the pan. After a while, 200.0000 g will be displayed on the screen. Place weight(s) of 200 g on the pan.</p>	 <p>PIONEER™</p> <p>200.0000 g</p> <p>Zero  Yes</p> <p>Print Unit No</p> <p>Function Mode Back</p> <p>Tare Menu-Cal Exit</p>
<p>6. Once the linearity calibration is completed successfully, CALdone will display. Press any button to return to the previous screen.</p>	 <p>PIONEER™</p> <p>CALdone</p> <p>Zero  Yes</p> <p>Print Unit No</p> <p>Function Mode Back</p> <p>Tare Menu-Cal Exit</p>
<p>7. Remove the weight(s), and the reading will be set to zero.</p>	 <p>PIONEER™</p> <p>* 0.0000 g</p> <p>Zero  Yes</p> <p>Print Unit No</p> <p>Function Mode Back</p> <p>Tare Menu-Cal Exit</p>

5.3 Balance Setup

Enter this sub-menu to customize the balance functionality.

Note: The factory default settings are shown below in **bold**.

5.3.1 Filter Level

Set the amount of signal filtering.		F I L T E R
Low = faster stabilization time with less stability. Medium = normal stabilization time with normal stability. High = slower stabilization time with more stability.		
Low Low	Med Medium	HIGH High

5.3.2 AZT (Auto Zero Tracking)

Set the automatic zero tracking functionality.		A Z T Auto Zero Tracking
Off = disabled. 0.5d = display maintains zero up to a drift of 0.5 graduation per second. 1d = display maintains zero up to a drift of 1 graduation per second. 3d = display maintains zero up to a drift of 3 graduations per second.		
0.5 d 0.5 d	1 d 1 d	3 d 3 d

5.3.3 Auto Tare

Set the automatic Tare.		
When Automatic Tare is set to On, the first object placed on the pan will be deemed as a container and tared automatically.		
Off = disabled. On = enabled.		
A.TArE Auto Tare	OFF Off	ON On

5.3.4 Graduations

Set the displayed readability of the balance.

1 Division = standard readability.

10 Divisions = readability is increased by a factor of 10.

For example, if the standard readability is 0.01 g, selecting 10 Divisions will result in a reading of 0.1 g.

 Graduation	 1 Division	 10 Division
----------------	----------------	-----------------

5.3.5 Date Format

Set the current date format.

Date Format:
 YYYY/MM/DD
MM/DD/YYYY
 DD/MM/YYYY

Date Format

MM/DD/YYYY

DD/MM/YYYY

YYYY/MM/DD

5.3.6 Date Setup

Set the current date in the desired date format.

To set the current date, press **No** or **Back** to increase or decrease the value.

Date

For example, if the current date is 22nd June, 2017,

MM/DD/YY: 06.22.17

DD/MM/YY: 22.06.17

YY/MM/DD: 17.06.22

5.3.7 Time Format

Set the time format.

24hr
12hr

Time Format

24hr

12hr

5.3.8 Time Setup

Set the current time in the desired time format.	t .mT
To set the current time, press No or Back to increase or decrease the value.	08.00.00

5.3.9 Brightness

Set the brightness of the display. Medium High Low	br ight Brightness	mT d Medium
	H IGH High	mT d Low

5.3.10 Auto Dim

Set whether the balance automatically turns off the display backlight of the display.	Aut. d .mT	
Off = disabled. 10 minutes = become dim if there is no motion in the display for 10 minutes. 20 minutes = become dim if there is no motion in the display for 20 minutes. 30 minutes = become dim if there is no motion in the display for 30 minutes.		
10 mT in 10 min	20 mT in 20 min	30 mT in 30 min

5.3.11 Approved Mode

Use this menu to set the Legal for Trade status. OFF = standard operation. ON = operation complies with Legal Metrology regulations.	L F t
---	-------

Note: When Approved Mode is set to On, the menu settings are affected as follows:

Calibration Menu:

- For InCal models, only Internal Calibration is available. All other functions are hidden.
- For ExCal models, the entire Calibration menu is hidden.

Balance Setup Menu:

- Filter Level is locked at the current setting.
- Auto Zero Tracking is limited to 0.5 Division and Off. The selected setting is locked.
- Auto Tare is locked at current setting.
- Graduations are forced to 1 Division and the menu item is hidden.

Communication Menu (Communication->Print Settings->Print Output):

- Stable Weight Only is locked On.
- Numeric Value Only is locked Off.

Communication Menu (Communication->Print Settings->Auto Print):

- Auto print mode selections are limited to Off, On Stability, and Interval. Continuous is not available.

Lockout Menu:

- Menu is hidden

Note: The security switch located at the rear of the balance must be in the locked position to set Approved Mode to On. The security switch must be in the unlocked position to set Approved Mode to Off. See Section 6.

5.4 Weighing Units

Enter this sub-menu to activate the desired units of measure.


PJX balances provide a choice of 9 units, which are all set On by default.

Note: Due to national laws, the balance may not include some of the units of measure listed below.

Unit	Display
Gram	g
Kilogram	kg
Milligram	mg
Carat	ct
Pound	lb
Ounce	oz
Ounce Troy	ozt
Pennyweight	dwt
Grain	GN

Changing Weighing Units

<p>1. Press and hold the Print / Unit button until the <i>Unit Menu</i> is displayed.</p>	
--	--

<p>2. The default unit is gram (g). To change the unit, press No to advance to the next unit.</p>	
<p>3. Press Yes to set the unit displayed to the weighing unit.</p>	

5.5 RS232 Interface Setup

<p>Enter this sub-menu to customize RS232 standard settings. Data may be output to either a printer or a PC.</p>	
--	--

5.5.1 Baud Rate

<p>Set the baud rate (bits per second).</p>	
<p>1200 = 1200 bps 2400 = 2400 bps 4800 = 4800 bps 9600 = 9600 bps 19200 = 19200 bps 38400 = 38400 bps</p>	

5.5.2 Transmission

Set the data bits, stop bit, and parity.

- 8-N-1 = 8 data bits, no parity, stop bit 1**
- 8-N-2 = 8 data bits, no parity, stop bit 2
- 7-E-1 = 7 data bits, even parity, stop bit 1
- 7-E-2 = 7 data bits, even parity, stop bit 2
- 7-N-1 = 7 data bits, no parity, stop bit 1
- 7-N-2 = 7 data bits, no parity, stop bit 2
- 7-O-1 = 7 data bits, odd parity, stop bit 1
- 7-O-2 = 7 data bits, odd parity, stop bit 2

<p>PARITY</p> <p>Parity</p>	<p>8-N-1</p> <p>8 data bits, stop bit 1</p>
------------------------------------	--

5.5.3 Handshake

Set the flow control method.

- NONE = no handshaking**
- XON-XOFF = XON/XOFF handshaking
- HARDWARE = hardware handshaking

<p>H.SHAKE</p> <p>Handshake</p>	<p>NONE</p> <p>None</p>
<p>XON-OFF</p> <p>Xon / Xoff</p>	<p>HARDW</p> <p>Hardware</p>

5.6 Print Settings

Enter this sub-menu to customize data transfer settings.

	<p>Print</p>
--	---------------------

5.6.1 Stable Only

Off = values are printed immediately regardless of stability.
 On = values are printed only when the stability criteria are met.

	<p>STABLE</p>
--	----------------------

5.6.2 Numeric Only

<p>Off = All selected results are printed. On = Only numeric data values are printed.</p>	<p>NUM</p>
---	------------

5.6.3 Single Header

<p>Off = Headers will be printed for every print requirement. On = Headers will be printed once a day.</p>	<p>5 in.HEAD</p>
--	------------------

5.6.4 Auto Print

<p>Enable or disable the functionality of auto print, and set the specific auto print mode.</p>	<p>A.Pr int</p>
---	-----------------

<p>1. Off = disabled</p>	<p>OFF</p>
---------------------------------	------------

<p>2. On Stability = printing occurs when the stability criteria are met.</p>	<p>ON.StAb</p>
--	----------------

<p>When On Stability is selected, set the conditions for printing. Load = Prints when the displayed load is stable.</p>	<p>LoAd</p>
---	-------------

<p>Load and Zero = Prints when the displayed load and zero reading is stable.</p>	<p>LoAd.ZEr</p>
---	-----------------

<p>3. Print Interval = printing occurs at the defined time interval. When Print Interval is selected, set the time interval using the numeric keypad. Note: Settings of 1 to 3600 seconds are available. Default is 0.</p>	<p>intEr</p>
---	--------------

<p>4. Continuous = printing occurs continuously.</p>	<p>Cont inu</p>
---	-----------------

5.6.5 Header

<p>On = the header is printed. Off = the header is not printed.</p>	<p>HEADER</p>
--	---------------

5.6.6 Date and Time

<p>On = the date and the time are printed. Off = neither the date nor the time is printed.</p>	<p>DATE</p>
---	-------------

5.6.7 Balance ID

<p>On = the balance ID is printed. Off = the balance ID is not printed.</p>	<p>BAL ID</p>
--	---------------

5.6.8 Balance Name

<p>On = the balance name is printed. Off = the balance name is not printed.</p>	<p>BAL.NAME</p>
--	-----------------

5.6.9 User Name

<p>On = the user name is printed. Off = the user name is not printed.</p>	<p>USR.NAME</p>
--	-----------------

5.6.10 Project Name

<p>On = the project name is printed. Off = the project name is not printed.</p>	<p>PRJ.NAME</p>
--	-----------------

5.6.11 Application Name

<p>On = the application name is printed. Off = the application name is not printed.</p>	<p>APP.NAME</p>
--	-----------------

5.6.12 Result

<p>On = the weighing result is printed. Off = the weighing result is not printed.</p>	<p>RESULT</p>
--	---------------

5.6.13 Gross

<p>On = the gross weight is printed. Off = the gross weight is not printed.</p>	<p>Gross</p>
---	--------------

5.6.14 Net

<p>On = the net weight is printed. Off = the net weight is not printed.</p>	<p>NET</p>
---	------------

5.6.15 Tare

<p>On = the tare weight is printed. Off = the tare weight is not printed.</p>	<p>TARE</p>
---	-------------

5.6.16 Line Feed

<p>1 Line = move the paper up one line after printing. 4 Lines = move the paper up four lines after printing.</p>	<p>FEED</p>
<p>1 LINE 1 Line</p>	<p>4 LINES 4 Lines</p>

5.7 GLP

<p>Enter this menu to set the Good Laboratory Practices (GLP).</p>	<p>GLP</p>
--	------------

5.7.1 Header

<p>Enables the printing of GLP headings. There are up to 3 headings available. Alphanumeric settings up to 16 characters are available for each header setting.</p>	<p>HEAdEr 1 Header 1</p>
<p>HEAdEr 2 Header 2</p>	<p>HEAdEr 3 Header 3</p>

5.7.2 Balance Name

Set the balance name.
Alphanumeric settings up to 16 characters are available for each Header setting.

BAL.NAM

5.7.3 User Name

Set the user name.
Alphanumeric settings up to 16 characters are available for each Header setting.

USR.NAM

5.7.4 Project Name

Set the user name.
Alphanumeric settings up to 16 characters are available for each Header setting.
The default is blank.

PRJ.NAM

5.8 Factory Reset

Use this sub-menu to reset the all menu settings to their Factory default settings.

Reset All = resets all menus to their factory default settings.
Exit = return to application main screen without resetting any menus.

RESET

5.9 Lockout

Use this sub-menu to lock / unlock certain menus.

Off = the menu is unlocked.
On = the menu is locked.

LOCK

6. LEGAL FOR TRADE (LFT)

When the balance is used in trade or a legally controlled application it must be set up, verified and sealed in accordance with local weights and measures regulations. It is the responsibility of the purchaser to ensure that **all pertinent legal requirements** are met.

6.1 Settings

Before verification and sealing, perform the following steps in order:

1. Verify that the menu settings meet the local weights and measures regulations.
2. Verify the units turned On meet the local weights and measures regulations.
3. Perform a calibration as explained in Section 5.2.
4. Set Approved Mode to On in the Balance Setup menu.
5. As shown in Figure 6-1, press the push button inside the hole.

Note: When Approved Mode is set to On, external calibration can't be performed.

6.2 Verification

A weights and measures official or authorized service agent must perform the verification procedure.

6.3 Sealing

After the Balance has been verified, it must be sealed to prevent undetected access to the legally controlled settings. Before sealing the device, ensure the Approved Mode setting in the Balance Setup menu has been set to ON.

- If using a paper seal, place seals over the security switch and the bottom housing as shown.
- If using a wire seal, pass the sealing wire through the holes in the security screw and the bottom housing as shown.

Figure 6-1. Sealing

7. PRINTING

7.1 Connecting, configuring and Testing the Printer / Computer Interface

Use the built-in RS-232 Port to connect either to a computer or a printer. If connecting to a computer, use HyperTerminal or similar software like SPDC described below.

(Find HyperTerminal under **Accessories/Communications** in Windows XP.)

Connect to the computer with a standard (straight-through) serial cable.


Choose **New Connection**, "connect using" COM1 (or available COM port).

Select **Baud=9600; Parity=8 None; Stop=1; Handshaking=None**. Click **OK**.

Choose Properties/Settings, then ASCII Setup. Check boxes as illustrated:

(Send line ends...; Echo typed characters...; Wrap lines...)


Use RS232 Interface Commands (Section 9.6.1) to control the balance from a PC.

SPDC Software

The Serial Port Data Collection / SPDC software is provided by Ohaus and can be used on operating systems that do not have the HyperTerminal software mentioned above. SPDC software can preliminarily collect and transfer the data to Microsoft files (such as Excel, Word, etc.).

Choose the export file type and export file path and then press Run as shown below.

System Requirements

- PC running Windows 98®, Windows 98SE®, Windows ME®, Windows 2000®, Windows XP®, Windows 7® or Windows 8® (32-bit).

Note: The latest SPDC software support English and Chinese language and can be downloaded from the Ohaus' website. For more information, refer to the *SPDC Data Collection Instruction Manual*.

7.2 Output format

The Result Data, and G/N/T data, is output in the following format.

Field:	Label ¹	Space ²	Weight ³	Space ²	Unit ⁴	Space	Stability ⁵	Space	G/N ⁶	Space	Term. Characters ⁷
Length:		1	11	1	5	1	≤ 1	≤ 1	≤ 3	0	≤ 8

Note:

1. The length of the label field is not fixed.
2. Each field is followed by a single delimiting space (ASCII 32).
3. The Weight field is 11 right justified characters. If the value is negative, the '-' character is located at the immediate left of the most significant digit.
4. The Unit field contains the unit of measure abbreviation up to 5 characters, right justified.
5. The Stability field contains the “?” character if the weight reading is not stable. The Stability field and the following Space field are omitted if the weight reading is stable.
6. The G/N field contains the net or gross indication. For net weights, the field contains “N”. For gross weights, the field contains “G”.
7. The Termination Characters field contains CRLF, Four CRLF or Form Feed (ASCII 12), depending on the LINE FEED menu setting.
8. When Numeric Only is set On, only the Weight Field is printed, left-aligned.

7.3 Printout Examples

Examples for each Application are shown with all items turned **ON** in the **Print** menu. The default values for **Header** lines 1-3 are also shown.

WEIGHING

```

Header 1
Header 2
Header 3
07/21/2017 17:41:44
Balance ID: B234567890
Balance Name: PJX8201/E
User Name: OHAUS
Project Name: PJX
Weigh
1365.1 g N
Gross: 3919.9 g G
Net: 1365.1 g NET
Tare: 2554.8 g T

Signature: _____
Verified By: _____
    
```

PART COUNTING

```

Header 1
Header 2
Header 3
07/21/2017 17:41:44
Balance ID: B234567890
Balance Name: PJX8201/E
User Name: OHAUS
Project Name: PJX
Count
Quantity: 201 PCS N
Gross: 300.1 g G
Net: 200.3 g N
Tare: 99.8 g T
APW: 1.00 g
Sample Size: 100 PCS

Signature: _____
Verified By: _____
    
```

PERCENT WEIGHING

Header 1
Header 2
Header 3
07/21/2017 17:41:44
Balance ID: B234567890
Balance Name: PJX8201/E
User Name: OHAUS
Project Name: PJX
Percent
Percentage: 20,025 % N
Gross: 300.1 g G
Net: 200.3 g N
Tare: 99.8 g T
Reference weight: 1000.1 g

Signature: _____
Verified By: _____

DENSITY DETERMINATION

Header 1
Header 2
Header 3
07/21/2017 17:41:44
Balance ID: B234567890
Balance Name: PJX8201/E
User Name: OHAUS
Project Name: PJX
Density
Density: 1.9943 g/cm3
Gross: 100.1 g G
Net: 100.1 g N
Tare: 0.0 g T
Weight in air: 200.1 g
Weight in water: 100.0 g

Signature: _____
Verified By: _____

INTERNAL CALIBRATION

-OHAUS-
07/26/2017 05:16:53
Balance ID:
Balance Name: PR2202
User Name:
Project Name:
---Internal Calibration---
Calibration is done.
Difference weight: 0.00 g

Signature: _____
Verified By: _____

SPAN CALIBRATION

-OHAUS-
07/26/2017 05:16:37
Balance ID:
Balance Name: PR2202
User Name:
Project Name:
---Span Calibration---
Calibration is done.
Reference weight: 2000.00 g
Actual weight: 2000.22 g
Difference weight: 0.22 g
Weight ID: _____

Signature: _____
Verified By: _____

LINEARITY CALIBRATION

-OHAUS-
07/26/2017 05:16:11
Balance ID:
Balance Name: PR2202
User Name:
Project Name:
---Linear Calibration---
Calibration is done.

Signature: _____
Verified By: _____

8. MAINTENANCE

8.1 Calibration

Periodically verify calibration by placing an accurate weight on the balance and viewing the result. If calibration is required, refer to section 5.2 for instructions.

8.2 Cleaning

WARNING: Disconnect the balance from the power supply before cleaning. Make sure that no liquid enters the interior of the balance.

Clean the Balance at regular intervals.

Housing surfaces may be cleaned with a lint-free cloth slightly dampened with water or a mild cleaning agent.

Glass surfaces may be cleaned with a commercial glass cleaner.

Attention: Do not use solvents, harsh chemicals, ammonia or abrasive cleaning agents.

8.3 Troubleshooting

TABLE 8-1 TROUBLESHOOTING

Symptom / Display	Possible Cause
<i>Err 8.1</i>	Weight reading exceeds Power On Zero limit.
<i>Err 8.2</i>	Weight reading below Power On Zero limit.
<i>Err 8.3</i>	Weight reading exceeds overload limit.
<i>Err 8.4</i>	Weight reading below underload limit.
<i>REF Err</i>	Reference Weight too small. The weight on the pan is too small to define a valid reference weight. Add additional samples.
<i>LowRef</i>	A warning that the reference weight is below the recommended value for accurate parts counting. Add additional samples if possible otherwise counting continues with diminished accuracy.
-----	Busy (tare, zero, and printing).
--RQ--	Function not executed.
<i>Error</i>	Timeout
<i>Idnr.Err</i>	IDNR Error

8.4 Service Information

If the troubleshooting section does not resolve your problem, contact an Authorized Ohaus Service Agent. Please visit our website www.ohaus.com to locate the Ohaus office nearest you.

9. TECHNICAL DATA

9.1 Specifications

Ambient conditions

- Indoor use only
- Altitude: Up to 2000 m
- Specified Temperature range: 10°C to 30°C
- Humidity: maximum relative humidity 80% for temperatures up to 30°C, decreasing linearly to 50% relative humidity at 40°C
- Operability is assured at ambient temperatures between 5°C and 40°C
- Mains supply voltage fluctuations: up to $\pm 10\%$ of the nominal voltage
- Installation category II
- Pollution degree: 2
- Supply voltage: 12V–0.5A

Materials

- Bottom Housing: Die-cast Aluminum, Painted + Plastic (HIPS)
- Top Housing: Plastic (HIPS)
- Weighing Platforms: Stainless steel
- Draft Shield: Glass, plastic (HIPS)
- Feet: Plastic (ABS)

TABLE 9-1. SPECIFICATIONS

Model	PJX203 PJX203/E	PJX803 PJX803/E	PJX1203 PJX1203/E	PJX1603 PJX1603/E	PJX2503 PJX2503/E
Capacity (ct)	200	800	1200	1600	2500
Readability d (ct)	0.001	0.001	0.001	0.005	0.005
Capacity (g)	40	160	240	320	500
Readability d (g)	0.0001	0.0001	0.0001	0.001	0.001
Repeatability (STDEV) (g)	0.0001	0.0001	0.0001	0.001	0.001
Linearity (g)	0.0002	0.0002	0.0002	0.002	0.002
Stabilization Time Typical (s)	3	3	3	3	3
Sensitivity Temperature Drift (PPM/K)	±3	±3	±3	±3	±3
Units	Milligram, Gram, Kilogram, Ounce, Pound, Carat, Pennyweight, Ounce Troy				
Applications	Basic Weighing; Parts Counting; Percent Weighing; Density Determination				
Platform Size (diameter)	3.5 in / 9 cm	3.5 in / 9 cm	3.5 in / 9 cm	4.7 in / 12 cm	4.7 in / 12 cm
Span Calibration Points (g)	20, 40	80, 160	100, 200	150, 300	250, 500
Linearity Calibration Points (g)	10, 20, 30, 40	40, 80, 120, 160	50, 100, 150, 200	100, 200, 300	200, 300, 400, 500
Tare Range	To capacity by subtraction				
Power Supply	Power input: 100-240V ~ 200mA 50-60Hz 12-18VA Power output: 12 VDC 0.5A				
Assembled Dimensions (W x D x H)	209 x 321 x 274 mm 8.2 x 12.6 x 10.8 inch				
Communication	RS232, MiniDin				
Operating Temperature Range	10°C to 30°C				
Storage Temperature Range	Humidity: maximum relative humidity 80% for temperatures up to 30°C, decreasing linearly to 50% relative humidity at 40°C				
Storage Conditions	-10°C to 60°C, humidity 10% to 90%, without condensation				
Net Weight	10 lb / 4.5 kg				
Shipping Weight	15.4lb / 7 kg				
Shipping Dimensions (W x D x H)	507 x 387 x 531 mm 20 x 15 x 21 inch				

*SRP refers to the standard deviation for n replicate weighings (n ≥ 10).

TABLE 9-2. SPECIFICATIONS (continued)

Model	PJX822 PJX822/E	PJX1602 PJX1602/E	PJX2202 PJX2202/E	PJX3202 PJX3202/E	PJX4202 PJX4202/E
Capacity (ct)	4100	8000	11000	16000	21000
Readability d (ct)	0.05	0.05	0.05	0.05	0.05
Capacity (g)	820	1600	2200	3200	4200
Readability d (g)	0.01	0.01	0.01	0.01	0.01
Repeatability (STDEV) (g)	0.01	0.01	0.01	0.01	0.01
Linearity (g)	0.02	0.02	0.02	0.02	0.02
Stabilization Time Typical (s)	2	2	2	2	2
Sensitivity Temperature Drift (PPM/K)	±6	±6	±6	±3	±3
Units	Milligram, Gram, Kilogram, Ounce, Pound, Carat, Pennyweight, Ounce Troy				
Applications	Basic Weighing; Parts Counting; Percent Weighing; Density Determination				
Platform Size (diameter)	7.1 in / 18 cm	7.1 in / 18 cm	7.1 in / 18 cm	7.1 in / 18 cm	7.1 in / 18 cm
Span Calibration Points (g)	400, 800	750, 1500	1000, 2000	1500, 3000	2000, 4000
Linearity Calibration Points (g)	200, 400, 600, 800	500, 1000, 1500	500, 1000, 1500, 2000	1000, 2000, 3000	1000, 2000, 3000, 4000
Tare Range	To capacity by subtraction				
Power Supply	Power input: 100-240V ~ 200mA 50-60Hz 12-18VA Power output: 12 VDC 0.5A				
Assembled Dimensions (W x D x H)	209 x 321 x 98 mm 8.2 x 12.6 x 3.9 inch				
Communication	RS232, MiniDin				
Operating Temperature Range	10°C to 30°C				
Storage Temperature Range	Humidity: maximum relative humidity 80% for temperatures up to 30°C, decreasing linearly to 50% relative humidity at 40°C				
Storage Conditions	-10°C to 60°C, humidity 10% to 90%, without condensation				
Net Weight	7.7 lb / 3.5 kg				
Shipping Weight	11 lb / 5 kg				
Shipping Dimensions (W x D x H)	550 x 385 x 291 mm 22 x 15 x 12 inch				

*SRP refers to the standard deviation for n replicate weighings (n ≥ 10).

TABLE 9-3. SPECIFICATIONS (continued)

Model	PJX5202 PJX5202/E	PJX2201 PJX2201/E	PJX4201 PJX4201/E	PJX6201 PJX6201/E
Capacity (ct)	26000	11000	21000	31000
Readability d (ct)	0.05	0.5	0.5	0.5
Capacity (g)	5200	2200	4200	6200
Readability d (g)	0.01	0.1	0.1	0.1
Repeatability (STDEV) (g)	0.01	0.1	0.1	0.1
Linearity (g)	0.02	0.2	0.2	0.2
Stabilization Time Typical (s)	2	1	1	1
Sensitivity Temperature Drift (PPM/K)	±3	±6	±3	±10
Units	Milligram, Gram, Kilogram, Ounce, Pound, Carat, Pennyweight, Ounce Troy			
Applications	Basic Weighing; Parts Counting; Percent Weighing; Density Determination			
Platform Size (diameter)	7.1 in / 18 cm	7.1 in / 18 cm	7.1 in / 18 cm	7.1 in / 18 cm
Span Calibration Points (g)	2500, 5000	1000, 2000	2000, 4000	3000, 6000
Linearity Calibration Points (g)	2000, 3000, 4000, 5000	500, 1000, 1500, 2000	1000, 2000, 3000, 4000	2000, 4000, 6000
Tare Range	To capacity by subtraction			
Power Supply	Power input: 100-240V ~ 200mA 50-60Hz 12-18VA Power output: 12 VDC 0.5A			
Assembled Dimensions (W x D x H)	209 x 321 x 98 mm 8.2 x 12.6 x 3.9 inch			
Communication	RS232, MiniDin			
Operating Temperature Range	10°C to 30°C			
Storage Temperature Range	Humidity: maximum relative humidity 80% for temperatures up to 30°C, decreasing linearly to 50% relative humidity at 40°C			
Storage Conditions	-10°C to 60°C, humidity 10% to 90%, without condensation			
Net Weight	7.7 lb / 3.5 kg			
Shipping Weight	11 lb / 5 kg			
Shipping Dimensions (W x D x H)	551 x 385 x 291 mm 22 x 15 x 12 inch			

*SRP refers to the standard deviation for n replicate weighings (n ≥ 10).

TABLE 9-4. SPECIFICATIONS (continued)

M model	PJX203M	PJX803M	PJX1203M	PJX1603M	PJX2503M
Capacity (ct)	200	800	1200	1600	2500
Readability d (ct)	0.001	0.001	0.001	0.01	0.01
Verification Interval e (ct)	0.01	0.01	0.01	0.1	0.1
Capacity (g)	40	160	240	320	500
Readability d (g)	0.0001	0.0001	0.0001	0.001	0.001
Verification Interval e (g)	0.001	0.001	0.001	0.01	0.01
Class	II	I	I	II	II
Repeatability (STDEV) (g)	0.0001	0.0001	0.0001	0.001	0.001
Linearity (g)	0.0002	0.0002	0.0002	0.002	0.002
Stabilization Time Typical (s)	3	3	3	3	3
Sensitivity Temperature Drift (PPM/K)	±3	±3	±3	±3	±3
Units	g, mg, ct				
Applications	Basic Weighing; Parts Counting; Percent Weighing; Density Determination				
Platform Size (diameter)	3.5 in / 9 cm			4.7 in / 12 cm	
Span Calibration Points (g)	20, 40	80, 160	100, 200	150, 300	250, 500
Linearity Calibration Points (g)	10, 20, 30, 40	40, 80, 120, 160	50, 100, 150, 200	100, 200, 300	200, 300, 400, 500
Tare Range	To capacity by subtraction				
Power Supply	Power input: 100-240V ~ 200mA 50-60Hz 12-18VA Power output: 12 VDC 0.5A				
Assembled Dimensions (W x D x H)	209 x 321 x 274 mm 8.2 x 12.6 x 10.8 inch				
Communication	RS232, MiniDin				
Operating Temperature Range	10°C to 30°C				
Storage Temperature Range	Humidity: maximum relative humidity 80% for temperatures up to 30°C, decreasing linearly to 50% relative humidity at 40°C				
Storage Conditions	-10°C to 60°C, humidity 10% to 90%, without condensation				
Net Weight	10 lb / 4.5 kg				
Shipping Weight	15.4lb / 7 kg				
Shipping Dimensions (W x D x H)	507 x 387 x 531 mm 20 x 15 x 21 inch				

*SRP refers to the standard deviation for n replicate weighings (n ≥ 10).

TABLE 9-5. SPECIFICATIONS (continued)

M model	PJX822M	PJX1602M	PJX2202M	PJX3202M	PJX4202M	PJX5202M	PJX4201M
Capacity (ct)	4100	8000	11000	16000	21000	26000	21000
Readability d (ct)	0.5	0.5	0.5	0.5	0.5	0.5	0.5
Verification Interval e (ct)	0.5	0.5	0.5	0.5	0.5	0.5	0.5
Capacity (g)	820	1600	2200	3200	4200	5200	4200
Readability d (g)	0.01	0.01	0.01	0.01	0.01	0.01	0.1
Verification Interval e (g)	0.1	0.1	0.1	0.1	0.1	0.1	0.1
Class	II	II	II	II	II	II	II
Repeatability (STDEV) (g)	0.01	0.01	0.01	0.01	0.01	0.01	0.1
Linearity (g)	0.02	0.02	0.02	0.02	0.02	0.02	0.2
Stabilization Time Typical (s)	2	2	2	2	2	2	1
Sensitivity Temperature Drift (PPM/K)	±6	±6	±6	±3	±3	±3	±3
Units	g, mg, ct						
Applications	Basic Weighing; Parts Counting; Percent Weighing; Density Determination						
Platform Size (diameter)	7.1 in / 18 cm						
Span Calibration Points (g)	400, 800	750, 1500	1000, 2000	1500, 3000	2000, 4000	2500, 5000	2000, 4000
Linearity Calibration Points (g)	200, 400, 600, 800	500, 1000, 1500	500, 1000, 1500, 2000	1000, 2000, 3000	1000, 2000, 3000, 4000	2000, 3000, 4000, 5000	1000, 2000, 3000, 4000
Tare Range	To capacity by subtraction						
Power Supply	Power input: 100-240V ~ 200mA 50-60Hz 12-18VA Power output: 12 VDC 0.5A						
Assembled Dimensions (W x D x H)	209 x 321 x 98 mm 8.2 x 12.6 x 3.9 inch						
Communication	RS232, MiniDin						
Operating Temperature Range	10°C to 30°C						
Storage Temperature Range	Humidity: maximum relative humidity 80% for temperatures up to 30°C, decreasing linearly to 50% relative humidity at 40°C						
Storage Conditions	-10°C-60°C, humidity 10% to 90%, without condensation						
Net Weight	7.7 lb / 3.5 kg						
Shipping Weight	11 lb / 5 kg						
Shipping Dimensions (W x D x H)	550 x 385 x 291 mm 22 x 15 x 12 inch						

*SRP refers to the standard deviation for n replicate weighings (n ≥ 10).

9.2 Drawings and Dimensions

Fully assembled dimensions


Figure 9-1 0.001 g / 0.0001 g model

Figure 9-2 0.1 g / 0.01 g model

9.3 Accessories

TABLE 9-6. ACCESSORIES

DESCRIPTION	PART NUMBER
Auxiliary Display AD7-MD	30472063
Security Device	80850043
Density Determination Kit	80253384
RS232 Cable (25 pin)	80500524
RS232 Cable (9 pin)	80500525
RS232 Interface	30268982
Dust Cover	30093334
In-use Cover	30372546
Printer SF40A	30064202 (EU); 30064203 (AM)
Power Adapter for Balance	46001724


9.4 Communication

9.4.1 Interface Commands

Commands listed in the following table will be acknowledged by the balance.

Command	Function
IP	Immediate Print of displayed weight (stable or unstable).
P	Print displayed weight (stable or unstable).
CP	Continuous Print.
SP	Print on Stability.
SLP	Set Auto Print to On Stability, allow non-zero displayed weight be printed.
SLZP	Set Auto Print to On Stability, allow both stable non-zero weight and stable zero reading printed.
xP	Set Auto Print to Interval Print, x = print interval (1-3600 sec), 0P disable the interval Print
0P	0P disable interval print, continuous print or print on stability
H	Enter Print Header Lines, the format is: H x "header string". Where x = line number 1 to 3, "header string" can be up to 24 alphanumeric characters. If no string in the command, "H x" will read the stored header x.
Z	Same as pressing Zero Key.
T	Same as pressing Tare Key.
xT***	Establish a preset Tare value in displayed unit. x = preset tare value. Sending 0T clears tare (if allowed).
PT	Prints Tare weight stored in memory.
PM	Print current application mode (weighing mode).
xM	Set current application mode to x. x depends on applications 1 – Weigh 2 – Count 3 – Percent 4 – Density
M	Scroll to the next enabled mode.
PU	Print Current weighing unit: g, kg, lb, oz, etc....
xU	Set balance to unit x: g, kg etc. 1 – g 2 – kg 3 – mg 4 – ct ...
U	Scroll to the next enabled unit.
ON	Brings out of Standby
OFF	Goes to Standby.
C	Begin Span Calibration
IC	Begin Internal Calibration.
AC	Abort Calibration.
PSN	Print Serial Number.
PV	Print Version: print name, software revision and LFT ON (if LFT is set ON).
x#	Set Counting APW (x) in grams. (must have APW stored)
P#	Print Counting application APW.
x%	Set Percent application reference weight (x) in grams. (must have reference weight stored)
P%	Print Percent application reference weight.
PTIME	Print current time.
PDATE	Print current date.
xTIME	Set Time x format: hh:mm:ss
xDATE	Set Date x format: mm/dd/yyyy
xS	0 = print unstable data, 1 = print stable only
xRL	Enable or disable OK response to non-print commands: x=0 to disable, x=1 to enable.
xT	Pre-tare the container weight (x) in grams.

9.4.2 RS232 (DB9) Pin Connections

Diagram	Type	Description
 <p>The diagram shows a 9-pin female DB9 connector with the following connections:</p> <ul style="list-style-type: none"> Pin 5: GND Pin 2: RxID (DATA IN) Pin 3: TxID (DATA OUT) Pin 4: RT_CTL (DATA OUT) Pin 9: GND Pin 8: CTS (HANDSHAKE IN) Pin 7: RTS (HANDSHAKE OUT) 	Interface type	Voltage interface conforming to EIA RS-232C/DIN 66020 (CCITT V24/V.28)
	Max. cable length	15 m
	Signal level	Output: +5 V ... +15 V (RL = 3 – 7kΩ) -5 V ... -15 V (RL = 3 - 7 kΩ) Input: +3 V ... +25 V -3 V ... -25 V
	Connector	Sub-D, 9-pole, female
	Operating mode	Full duplex
	Transmission mode	Bit-serial, asynchronous
	Transmission code	ASCII
	Baud rates	1200, 2400, 4800, 9600, 19200, 38400 (firmware selectable)
	Bits/parity	7-bit/even, 7-bit/odd, 7-bit/none, 8-bit/none (firmware selectable)
	Stop bits	Stop bit 1, 2
	Handshake	None, XON/XOFF, RTS/CTS (selectable)
	End-of-line	Not selectable

10. SOFTWARE UPDATES

Ohaus is continuously improving its balance software. To obtain the latest release, please contact your Authorized Ohaus Dealer or Ohaus Corporation.

11. COMPLIANCE

Compliance to the following standards is indicated by the corresponding mark on the product.


Mark	Standard
	This product complies with the EU Directives 2011/65/EU (RoHS), 2014/30/EU (EMC), 2014/35/EU (LVD) and 2014/31/EU (NAWI). The EU Declaration of Conformity is available online at www.ohaus.com/ce .
	EN 61326-1, AS/NZS 61010-1
	CAN/CSA-C22.2 No. 61010-1 UL Std. No. 61010-1

	This product complies with the EU Directive 2012/19/EU (WEEE). Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. For disposal instructions in Europe, refer to www.ohaus.com/wEEE .
--	---

Important notice for verified weighing instruments in the EU

When the instrument is used in trade or a legally controlled application it must be set up, verified and sealed in accordance with local weights and measures regulations. It is the responsibility of the purchaser to ensure that all pertinent legal requirements are met.

Weighing Instruments verified at the place of manufacture bear the following supplementary metrology marking on the descriptive plate.

Weighing Instruments to be verified in two stages have no supplementary metrology marking on the descriptive plate. The second stage of conformity assessment must be carried out by the applicable weights and measures authorities.

If national regulations limit the validity period of the verification, the user of the weighing instrument must strictly observe the re-verification period and inform the weights and measures authorities.

As verification requirements vary by jurisdiction, the purchaser should contact their local weights and measures office if they are not familiar with the requirements.

FCC Note

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Industry Canada Note

This Class A digital apparatus complies with Canadian ICES-003.

ISO 9001 Registration

In 1994, OHAUS Corporation, USA, was awarded a certificate of registration to ISO 9001 by Bureau Veritas Quality International (BVQI), confirming that the OHAUS quality management system is compliant with the ISO 9001 standard's requirements. On June 21, 2012, OHAUS Corporation, USA, was re-registered to the ISO 9001:2008 standard.

LIMITED WARRANTY

Ohaus products are warranted against defects in materials and workmanship from the date of delivery through the duration of the warranty period. During the warranty period Ohaus will repair, or, at its option, replace any component(s) that proves to be defective at no charge, provided that the product is returned, freight prepaid, to Ohaus.

This warranty does not apply if the product has been damaged by accident or misuse, exposed to radioactive or corrosive materials, has foreign material penetrating to the inside of the product, or as a result of service or modification by other than Ohaus. In lieu of a properly returned warranty registration card, the warranty period shall begin on the date of shipment to the authorized dealer. No other express or implied warranty is given by Ohaus Corporation. Ohaus Corporation shall not be liable for any consequential damages.

As warranty legislation differs from state to state and country to country, please contact Ohaus or your local Ohaus dealer for further details.

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Calibration	Setup	Unit	RS232	Print	GLP	Reset	Lock
InCal	Filter Level	Gram	Baud Rate	Stable Only	Header 1	Reset All	Calibration
Cal Adjust	AZT	Kilogram	Parity	Numeric Only	Header 2		Setup
Span Cal	Auto Tare	Milligram	Handshake	Single Header	Header 3		RS232
Linearity Cal	Graduations	Carat		Auto Print	Balance Name		Print
	Date Format	Pound		Header	User Name		GLP
	Date Setting	Ounce		Date and Time	Project Name		Reset
	Time Format	Ounce Troy		Balance ID			
	Time Setting	Pennyweight		Balance Name			
	Brightness	Grain		User Name			
	Auto Dim			Project Name			
	LFT			Application Name			
			Result				
			Gross Weight				
			Net Weight				
			Tare Weight				
			Line Feed				

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PESAGE

```
Header 1
Header 2
Header 3
07/21/2017 17:41:44
Balance ID: B234567890
Balance Name: PJX8201/E
User Name: OHAUS
Project Name: PJX
Weigh
1365.1 g N
Gross: 3919.9 g G
Net: 1365.1 g NET
Tare: 2554.8 g T

Signature: _____
Verified By: _____
```

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```
Header 1
Header 2
Header 3
07/21/2017 17:41:44
Balance ID: B234567890
Balance Name: PJX8201/E
User Name: OHAUS
Project Name: PJX
Count
Quantity: 201 PCS N
Gross: 300.1 g G
Net: 200.3 g N
Tare: 99.8 g T
APW: 1.00 g
Sample Size: 100 PCS

Signature: _____
Verified By: _____
```

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Header 1
 Header 2
 Header 3
 07/21/2017 17:41:44
 Balance ID: B234567890
 Balance Name: PJX8201/E
 User Name: OHAUS
 Project Name: PJX
 Percent
 Percentage: 20,025 % N
 Gross: 300.1 g G
 Net: 200.3 g N
 Tare: 99.8 g T
 Reference weight: 1000.1 g

Signature: _____
 Verified By: _____

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Header 1
 Header 2
 Header 3
 07/21/2017 17:41:44
 Balance ID: B234567890
 Balance Name: PJX8201/E
 User Name: OHAUS
 Project Name: PJX
 Density
 Density: 1.9943 g/cm3
 Gross: 100.1 g G
 Net: 100.1 g N
 Tare: 0.0 g T
 Weight in air: 200.1 g
 Weight in water: 100.0 g

Signature: _____
 Verified By: _____

Calibrage interne

-OHAUS-
 07/26/2017 05:16:53
 Balance ID:
 Balance Name: PR2202
 User Name:
 Project Name:
 ---Internal Calibration---
 Calibration is done.
 Difference weight: 0.00 g

Signature: _____
 Verified By: _____

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-OHAUS-
 07/26/2017 05:16:37
 Balance ID:
 Balance Name: PR2202
 User Name:
 Project Name:
 ---Span Calibration---
 Calibration is done.
 Reference weight: 2000.00 g
 Actual weight: 2000.22 g
 Difference weight: 0.22 g
 Weight ID: _____

Signature: _____
 Verified By: _____

Calibrage de linéarité

-OHAUS-
 07/26/2017 05:16:11
 Balance ID:
 Balance Name: PR2202
 User Name:
 Project Name:
 ---Linear Calibration---
 Calibration is done.

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P/N 30372567 F © 2018 Ohaus Corporation, all rights reserved / todos los derechos reservados / tous droits réservés / alle Rechte vorbehalten / tutti i diritti riservati